



DrugSense[®]

PRODUCT INFO PACK



Andatech Corporation Pty Ltd

1300 800 200

andatech@andatech.com.au

www.andatech.com.au

Updated May 2014



ACUDRUG SALIVA DRUG TEST

Model No: DAIS-01 Bulk Carton of 25 RRP: \$420



It is a very easy test – simply stick the swab in the mouth and when saturated, typically 30 seconds or so, remove it and screw it into the test tube. Results are typically ready in less than 3 minutes start to finish making it one of the fastest and easiest devices to use on the market today. Additionally we have personally had excellent success with it and can speak to its accuracy in detecting drug use.

The THC cutoff is 75 ng/ml which is quite low and is for the parent drug.

The AcuDrug was specifically designed for the Australian market as it is set to the AS 4760 cutoff levels. AS 4760 provides the recommended cutoff levels for saliva testing in Australia. This is the only device in Australia set to these levels with the lone exception of THC. Additionally, it is important to

note that no device anywhere in the world meets the AS 4760 suggested levels for THC.

The AcuDrug device is also CASA compliant (Civil Aviation Safety Authority).

Features

- Set to AS 4760 target figures for all drug groups except THC
- Tests for the Parent THC
- Error Free – Collection test and swab in 1 device.
- Less invasive than urine
- Test anywhere – anytime
- From collection to result in less than 5 minutes
- Easy, Fast and Accurate 1-step test
- Screens for the Delta 9 ($\Delta 9$) THC Component, not the metabolite.
- Test for 6 drugs: THC / COC / AMP / METH / OPI / BZO

Test Drug	Cut-off level ($\mu\text{g/L}$)
AMP - Amphetamine	50 ng/ml
BZO - Benzodiazepines	50 ng/ml
COC - Cocaine Metabolites	50 ng/ml
METH - Methamphetamine	50 ng/ml
OPI - Opiates	50 ng/ml
THC - Cannabis Metabolites	75 ng/ml



Andatech Corporation Pty Ltd

1300 800 200

andatech@andatech.com.au

www.andatech.com.au

DrugSense®

PRODUCT INFO PACK

QUICK PROFILE DOA SPLIT CUP URINE DRUG TEST

Model No: DOA6+2+K2 Bulk Carton of 30 RRP: \$555



The Quick Profile DOA Split Cup Urine Drug Test is an excellent drug testing kit for the qualitative assessment of drugs and their metabolites in human urine. The DOA Split Cup Urine Drug Test is also made to Australian Standards 4308:2008 and is ideal for drug testing in the workplace.

The DOA Split Cup is a rapid diagnostic test for drugs in the urine. The sample is collected in the cup, which separates the specimen into two distinct chambers for testing.

One portion of the specimen is used for the initial screen, and the other portion is reserved in the event that there is a presumptive positive and confirmation is needed through

further testing at a reference laboratory.

The DOA Split Cup tests for Amphetamines and Methamphetamines (including ecstasy), Benzodiazepines, Cocaine, Opiates, Cannabinoids (THC) and Synthetic Cannabis (K2) - pH and Creatinine.

The Quick Profile DOA Split Cup Urine Drug Test is easy to use and provides accurate, on-site test results within 5 minutes.

Test Drug	Cut-off level (µg/L)
AMP - Amphetamine	300 ng/ml
BZO - Benzodiazepines	200 ng/ml
COC - Cocaine Metabolites	300 ng/ml
METH - Methamphetamine	300 ng/ml
OPI - Opiates	300 ng/ml
THC - Cannabis Metabolites	50 ng/ml
Synthetic Cannabis (K2)	50 ng/ml



ORA-CHECK 5 SALIVA DRUG TEST

Model No: DMNS-05PLUS Bulk Carton of 25 RRP: \$412.50



The Ora-Check 5 Saliva Drug Test is made to Australian Standards AS4670 requirements and tests for the 5 drug groups:

- Cannabinoids THC, Marijuana
- Opiates, Heroin, Morphine
- Cocaine
- Methamphetamine, Party drugs
- Amphetamines, stimulants

The Ora-Check 5 Saliva Drug Test is very easy to use with just a simple 3-step process. Results are processed within a few minutes using easy to read lines.

The testing is non-invasive and collection is user friendly. Test your employees anywhere and anytime with this compact, easy to use drug testing kit.

Test Drug	Cut-off level (µg/L)	Detection Time
AMP - Amphetamine	50 ng/ml	10 minutes - 72 hours
BZO - Benzodiazepines	20 ng/ml	10 minutes - 72 hours
COC - Cocaine Metabolites	20 ng/ml	10 minutes - 24 hours
METH - Methamphetamine	50 ng/ml	10 minutes - 72 hours
OPI - Opiates	40 ng/ml	1 hour - 72 hours
THC - Cannabis Metabolites	50 ng/ml	1 hour - 24 hours

K2 QUICK VIEW™ SYNTHETIC MARIJUANA URINE DRUG TEST

Model No: DMNS-QV-K2 Bulk Carton of 25 RRP: \$187.50



The K2 Quick View™ Synthetic Marijuana Urine Drug Test is CE registered and uses a urine cassette test with pipette. The drug testing kit needs just 3 drops of sample in the cassette in order to conduct a test.

The K2 Marijuana urine drug test is an immunochromatography-based one-step in-vitro test. It is designed for quantitative determination of synthetic cannabis JWH-18 and JWH-73's major metabolites in human urine specimens at cut off levels of 50ng/ml.

Synthetic cannabinoids are chemical compounds that mimic the effect of THC, the principle active ingredient of cannabis. Like THC, they bind to cannabinoid receptors in the brain and were initially developed as therapeutic agents for pain treatment. These psychoactive research chemicals

are frequently being used in hermal mixtures and sold as 'fake weed' or 'synthetic marijuana'. JWH-018 and JWH-073 are the two most common synthetic cannabinoid chemicals found in a variety of herbal smoking blends. Contrary to common belief, legislation is in place in Australia prohibiting the use of these substances. See the Downloads tab for further information on legislation and the K2 drug testing kit.

Test Drug	Cut-off level (µg/L)
JWH-018 pentanoic acid	50 ng/ml
JWH-073 butanoic acid	50 ng/ml

PLEASE NOTE: The K2 Quick View™ Synthetic Marijuana Urine Drug Test provides only a preliminary analytical test result. A more specific chemical method must be used to obtain a confirmed analytical result. Gas chromatography/ mass spectrometry (GC/MS) has been established as the preferred confirmatory method. Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

ORAL-VIEW™ 5+ SALIVA DRUG TEST

Model No: DMNS-05 Bulk Carton of 25 RRP: \$262.50



The Oral-View™ 5+ Saliva Drug Test is the best party drug testing kit. It's a one-step rapid qualitative immunoassay for screening potential abuse of one or more drugs in our saliva. The Oral-View 5+ tests for:

- Ice - Methamphetamines - MET
- Ecstasy - MDMA / XTC
- Speed - Amphetamines - AMP
- Marijuana - THC
- Heroin, Morphine - Opiates - OPI
- Cocaine - Coc
- Sedatives - Benzodiazepines - BZO

With new, lower cut-off levels, a pipette for direct oral samples, and made in the USA to Australian Standard AS4670 target requirements, the Oral-View™ 5+ Saliva Drug Test is an easy and non-invasive way to test employees or individuals for drugs.

Test Drug	Cut-off level (µg/L)
AMP - Amphetamine	50 ng/ml
COC - Cocaine Metabolites	20 ng/ml
METH - Methamphetamine	50 ng/ml
OPI - Opiates	40 ng/ml
THC - Cannabis Metabolites	12 ng/ml